



# CT Patient lines

CT patient lines are used to ensure patient safety when using contrast media syringes

Item no	Description	Packaging Box/Carton
CTPL020RV-H	CT Patient line 20cm, 350 psi	100/600
CTPL040RV-H	CT Patient line 40cm, 350 psi	100/600
CTPL060RV-H	CT Patient line 60cm, 350 psi	100/600
CTPL075RV-H	CT Patient line 75cm, 350 psi	100/600
CTPL120RV-H	CT Patient line 120cm, 350 psi	50/300
CTPL020RV-HR	CT Patient line 20cm, 350 psi, RA	100/600
CTPL040RV-HR	CT Patient line 40cm, 350 psi, RA	100/600
CTPL060RV-HR	CT Patient line 60cm, 350 psi, RA	100/600
CTPL075RV-HR	CT Patient line 75cm, 350 psi, RA	100/600
CTPL120RV-HR	CT Patient line 120cm, 350 psi, RA	50/300




Item no	Description	Packaging Box/Carton
CTPL140RV-HY	CT Patient line 140cm, 350 psi with integrated y-connector	100/600

# CT Patient lines

CT patient lines are used to ensure patient safety when using contrast media syringes

Item no	Description	Packaging Box/Carton
CTPL020RV2-H	CT Patient line 20cm, 350 psi, 2 valves	100/600
CTPL040RV2-H	CT Patient line 40cm, 350 psi, 2 valves	100/600
CTPL060RV2-H	CT Patient line 60cm, 350 psi, 2 valves	100/600
CTPL075RV2-H	CT Patient line 75cm, 350 psi, 2 valves	100/600
CTPL120RV2-H	CT Patient line 120cm, 350 psi, 2 valves	50/300
CTPL020RV2-HR	CT Patient line 20cm, 350 psi, 2 valves, RA	100/600
CTPL040RV2-HR	CT Patient line 40cm, 350 psi, 2 valves, RA	100/600
CTPL060RV2-HR	CT Patient line 60cm, 350 psi, 2 valves, RA	100/600
CTPL075RV2-HR	CT Patient line 75cm, 350 psi, 2 valves, RA	100/600
CTPL120RV2-HR	CT Patient line 120cm, 350 psi, 2 valves, RA	50/300



## Features



- Saves time
- Female Luer-Lock fits most CT injectors
- During the injection phase the contrast media is injected in the usual way
- The transparent tubing helps finding bubbles
- Male Luer-Lock at distal end
- Non-return valve prevents contamination
- The patient line is exchanged for each patient.

## Advantages

- Saves contrast media
- Ensures patient safety – prevents microbiological contamination
- Fits most CT injectors
- Saves time for preparations

# CT Filling Sets

CT Filling Sets are used to ensure safe and easy refilling of CT syringes



Item no	Description	Packaging Box/Carton
CTFS120DRV	CT Filling Set 120cm, 350 psi	50/300
CTFS150DRV	CT Filling Set 150cm, 300 psi with drip chamber	50/300

## Features

- During the refill phase the contrast media is pulled into the syringe
- Saves time
- The valve allows for refilling the syringe without disconnecting the patient line
- Female Luer-Lock fits most CT injectors
- During the injection phase the contrast media is injected in the usual way
- The transparent tubing helps finding bubbles
- The spike has an air filter to prevent contamination

## Advantages

- Saves contrast media
- Prevents microbiological contamination
- Fits most CT injectors
- Saves time for preparations